EXPAIN[®] ACTIVATE KNEE

USER MANUAL

- › English
- › Norsk
- › Dansk
- › Svensk
- › Suomi

EXPAIN[®]

Contents

- 1 Components
- 2 Introduction
- 3 Before use
- 4 Support function
- 5 TENS function
- 6 Heating function
- 7 Control with app
- 8 Batteries and charging
- 9 Maintenance and storage
- 10 Technical specifications
- 11 Safety information

Components

PARTS

- 1 Knee support
 - a. TENS electrodes
 - b. Heat area
 - c. Magnetic control unit connectors
- 2 Control unit
 - a. Magnetic connectors to knee strap
 - b. Micro USB charging connector
- 3 USB cable and charger
- 4 User manual

CONTROL

- 5 on / off
- 6 TENS intensity up
- 7 TENS intensity down
- 8 Heat function on / off
- 9 Indicator Bluetooth
- 10 Indicator TENS

Button	Function	Standby LED	Heat LED
On/Off (5)	On	Purple	
	Bluetooth connected	Blue	
Heating (8)	Turn the heating On/Off		Red flash, equivalent heat level
	Low battery	Red flash	
	Charger	Red	
	Charged	Purple	
TENS + (6)	Increase TENS intensity		Blue flash
TENS - (7)	Reduce TENS intensity		Blue flash
TENS + (6) hold 1 sec	Change TENS program up		Two blue flashes
TENS - (7) hold 1 sec	Change TENS program down		Two blue flashes

Introduction

Congratulations on your new **EXPAIN Activate Knee**. This is a unique, wireless and Bluetooth controlled unit with active heat, good support and TENS technology specifically designed for painful knees.

This innovative device is safe, without medication, easy to use, discreet and comfortable to wear - and above all it gives you ability to control your pain, so you can live an active life.

Before Use

EXPAIN Activate Knee is designed for personal use. Read the instructions and safety precautions before use.

You need to charge the control unit before use. Insert the USB cable into the charger and place the charger in the wall socket. Make sure the charger fits into the wall socket, then connect the cable to the controller. The indicator light on the controller is red when charging. The device is fully charged after approximately 2 hours (2000 mAh).

The support function

Support function on **EXPAIN Activate Knee** is analog and ready for use. To put on the knee support:

1. Make sure that the knee support is the correct way. Magnetic contacts should be outwards from the body and pointing up and forward on the thigh. The hole in the knee support shall be placed so the kneecap is visible through the hole.
2. Pull up the zipper so that the knee support wraps around the knee. Then insert the belt buckles.
3. The knee support can be adjusted by tightening or loosening the upper and lower mounting strap.

Mounting the control unit

The control unit attaches to the knee support with the four magnetic buttons. The controller is placed with control buttons out and upwards. Ensure that all four magnetic buttons are fixed before the product is turned on.

TENS function

Make sure the skin under knee strap is clean. Attach knee strap and controller as indicated. If TENS treatment is not strong enough, or it is perceived as unpleasant or pungent, it is recommended to moisten the electrode area with water. Take off controller before the electrode areas on the Knee support are lightly moistened. The unit is powered on by pressing the On / Off button on the main unit. TENS intensity is zero at startup. To turn on TENS, press + button on the controller. Adjust intensity to the desired level. For best effect set intensity so that you feel it, but not so high that it hurts. If intensity is painful, turn down the intensity immediately with the button, or detach the control unit from the body. Change TENS program with the same buttons, press and hold for about 1 second.

Heating function

The knee support and its insulating properties provides passive heating with immediate effect when it is put on. Put on the knee support as described in section "Support function." The knee support also has active heating elements embedded.

The Heat function is activated by pressing heat On / Off button on the controller. Turn the heat off again by pressing heat On / Off button until there are no red flashes on the heat button.

The Indicator light on the heat button indicates heating level and heat on or off with red flashes. Immediately turn off the heat function or take off knee strap if it becomes uncomfortably hot.

Control with EXPAIN app

You can get more out of **EXPAIN Activate Knee** by connecting it to your smartphone. You can follow when and how much you've spent using **EXPAIN Activate Knee**. The **EXPAIN app** lets you control TENS and heat, and access additional TENS programs.

EXPAIN Activate Knee is turned on and set in Bluetooth pairing mode when you turn on the unit. If an already paired phone is close, the device will automatically connect to that phone. It is not possible to connect the device to more phones simultaneously. For someone else to pair, the device must be close to the app and verify that the standby light is purple. The unit turns off by pressing the On / Off button.

The app is available for download from the Apple App Store or Google Play Store. Search **EXPAIN** to find the app. Follow the instructions in the app.

The app also has a built in pain logging function which lets you compare your use of **EXPAIN Activate Knee** with your pain. This way you can see relations with your pain development. For best insight into your **EXPAIN Activate Knee** use and improvement of pain, it is recommended to enter pain registrations regularly.

Battery and charging

EXPAIN Activate Knee is charged with supplied or equivalent USB cable.

The device has a battery life of approximately 45 minutes using active heat and TENS. Without the use of active heating the battery capacity equals approximately 50 TENS treatments. The device alerts low battery with red flashes on the standby button. Take unit off the body and recharge it fully before it is used again. Do not charge the device when it is on the body. Indicator light is red when charging.

When the device is fully charged the indicator light is purple, and it's ready for use again. The main unit is normally charged from empty to full in about 2 hours.

Maintenance and Storage

Clean the device with a dry or slightly damp cloth.

The control unit is safest in the accompanying box. Bring the USB charging cable in the box so you always have this available.

The knee support can be machine washed, at 40°C or by hand washing. Remember to take off the control unit before the knee support is washed.

Ensure that the control unit does not get wet. It is recommended that the control unit is not used in humid environments or if you plan to do activities that will involve sweating. The knee support can withstand moisture and water and can be used during activities that cause sweating.

The unit should be transported, stored

and operated at a temperature between 10 ° C and 40 ° C, at a relative humidity between 30% and 85% and with an atmospheric pressure between 500 and 1060 hPa.

Hang up the knee support to dry after use. Keep the control unit in its original box and away from children.

If the device is not working properly, stop all operation and contact support at www.expain.com

Technical specifications	
Number TENS channels	1
Active heating function	3W
Technology in use	TENS, active heat, support and passive heat (insulation)
Power	LiPo rechargeable battery (2000 mAh)
Control	Bluetooth 4.1 and Control buttons
Size (DxWxH)	6 cm x 15 cm x 23 cm
Weight	0,8 kg
Safety standards	IEC 62133, EN 60950
EMC	EN 301 489-1, EN 301 489-17, EN 300 328
Made in	China

Safety information

Applications

For painful knees.

Safety rules

Do not use this product if you have very sensitive skin.

Warning!

Do not use the device in other locations on the body than those specified in this manual.

Do not use TENS function if you have a pacemaker, implanted defibrillator or heart problems. Consult your doctor if you are unsure.

TENS activation by other activities can increase the risk of injuries. Do not use the TENS-function when you bathe, shower, sleep, run, control machinery or perform other activities that

can be disturbed by TENS activation.

*Using TENS function near measuring devices (such as ECG device, **EXPAIN Change Neck** or **EXPAIN Change Arm** and ECG alarms) may cause the device to stop working. Do not use TENS function near electrical measuring instruments.*

*The efficacy and safety using TENS function on children, pregnant women, or on the head has not been studied. Do not use the TENS function on children, or let children use **EXPAIN Activate Knee**. Do not use **EXPAIN Activate Knee** if you are pregnant, or suspect that you are pregnant, except if your therapist recommends it. Do not use **EXPAIN Activate Knee** on the head.*

If you have or suspect you have epilepsy or heart failure, TENS activation may have unexpected effects.

Always ask a doctor before using TENS activation if you have epilepsy or heart defects.

Expain AS
Karl Johans gate 13
0154 Oslo, Norway
www.expain.com

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